




Separation of Tocopherols on HALO[®] C30 based on GB (Chinese Standards)

Application Note 189-V



PEAK IDENTITIES:

1. δ -tocopherol
2. γ -tocopherol
3. β -tocopherol
4. α -tocopherol

Tocopherols are forms of vitamin E (fat-soluble) that have antioxidant properties in both the human body and in food. They are also used for cosmetics and many personal care products. Here, tocopherols are separated on a 250 mm 160 Å pore size HALO[®] C30 column using a GB (Chinese standard) method. Due to the shape selectivity of the C30 phase, separation of the four isomers is achieved.

TEST CONDITIONS:

Column: HALO 160 Å C30, 2.7 μ m,
4.6 x 250 mm

Part Number: 92114-930

Mobile Phase:

- A: Water
- B: Methanol

Isocratic: 95% B

Flow Rate: 0.9 mL/min

Initial Pressure: 240 bar

Temperature: 30 °C

Detection: UV 294 nm, PDA

Injection Volume: 20 μ L

Sample Solvent: Methanol

Response Time: 2.0 sec


Data Rate: 20 Hz

Flow Cell: 13 μ L

LC System: Agilent 1100

Data Courtesy of Beijing Institute for Drug Control

STRUCTURE:



Tocopherol	R1	R2
Alpha (α)	CH ₃	CH ₃
Beta (β)	CH ₃	H
Gamma (γ)	H	CH ₃
Delta (δ)	H	H

