




Effect of Temperature on the Separation of Proteins on HALO 400 Å C4

Application Note 103-PR



PEAK IDENTITIES:

1. Lysozyme (14.3 kDa)
2. Bovine serum albumin (66.4 kDa)
3. α -Chymotrypsinogen A (25.0 kDa)
4. Enolase (46.7 kDa)
5. Ovalbumin (44.0 kDa)

These separations demonstrate the effect of elevated temperatures on the efficiency of protein separations done under reversed-phase conditions on a HALO 400 Å C4, 3.4 μm , column. One observes larger and narrower peaks as the temperature increases. The HALO[®] C4 phase has been shown to be very stable even at these elevated temperatures.

TEST CONDITIONS:

Column: HALO 400 Å C4, 3.4 μm ,
2.1 x 100 mm

Part Number: 93412-614

Mobile Phase: 72/28 - A/B

A: 0.1% trifluoroacetic acid in water

B: 0.1% trifluoroacetic acid in acetonitrile

Gradient: 28% B to 58% B in 10 min

Gradient Delay Volume: ~250 μL

Flow Rate: 0.45 mL/min

Pressure: See chart

Temperature: See chart

Detection: UV 215 nm, PDA

Injection Volume: 2.0 μL

Sample Solvent: Mobile phase A

Response Time: 1.0 sec

Flow Cell: 2.0 μL micro cell

LC System: Agilent 1200 SL

