




## Separation of Four Small Proteins on HALO® 160 Å ES-C18, 5 µm vs. Totally Porous C18, 3.0 µm

Application Note 104-PR



### PEAK IDENTITIES:

1. Ribonuclease A (13.7 KDa)
2. Cytochrome c (12.4 KDa)
3. Lysozyme (14.3 KDa)
4. α-Lactalbumin (14.2 KDa)

These chromatograms show the separation of four low MW proteins on HALO 160 Å ES-C18, 5 µm column vs. a totally porous C18, 3.0 µm column. The separations are similar with the benefit of the HALO® 5 µm column having lower back pressure and similar resolution. The HALO® 5 µm ES-C18 phase is made with sterically hindered silanes during manufacture, enhancing the stability-even at temperatures up to 90 °C. The stability of the totally porous C18 column was not evaluated.

### TEST CONDITIONS:

#### Columns:

1) HALO 160 Å ES-C18, 5 µm, 4.6 x 150 mm

**Part Number:** 95124-702

2) 100 Å totally porous C18, 3.0 µm, 4.6 x 150 mm

**Mobile Phase:** 72/28 - A/B (start)

A: Water with 0.1% trifluoroacetic acid

B: Acetonitrile with 0.1% trifluoroacetic acid

**Gradient:** 28% B to 55% B in 5 min

**Flow Rate:** 1.5 mL/min

**Pressure:** 95 bar (HALO®)

170 bar (competitor)

**Temperature:** 60 °C

**Detection:** UV 280 nm, PDA

**Injection Volume:** 15 µL

**Sample Solvent:** Mobile phase A

**Response Time:** 0.1 sec

**Flow Cell:** 2.0 µL micro cell

**LC System:** Agilent 1200 SL

