




Effect of Silica Pore Size on Protein Separations

Application Note 130-PR



PEAK IDENTITIES:

1. Ribonuclease A (13.7 kDa)
2. Cytochrome C (12.4 kDa)
3. Lysozyme (14.3 kDa)
4. α -Lactalbumin (14.2 kDa)
5. Catalase (tetramer of ~60 kDa each)
6. Enolase (46.7 kDa)

Sharper, taller peaks are observed using the HALO 400 Å ES-C18 column because the larger pore size allows unrestricted diffusion for these biomolecules into and out of the porous shell. The half height peak widths above each protein peak are significantly smaller with the HALO 400 Å column despite the larger particle size of the packing material, emphasizing the importance of larger pores when separating proteins.

TEST CONDITIONS:

Columns:

- 1) HALO 160 Å ES-C18, 2.7 μm , 4.6 x 100 mm
Part Number: 92124-602
- 2) HALO 400 Å ES-C18, 3.4 μm , 4.6 x 100 mm
Part Number: 93414-602

Mobile Phase:

- A: 0.1% trifluoroacetic acid in water
B: 0.1% trifluoroacetic acid in acetonitrile

Gradient: 23% B to 50% B in 15 min

Flow Rate: 1.5 mL/min

Initial Pressure: See chart

Temperature: 60 °C

Detection: UV 215 nm, VWD

Injection Volume: 5.0 μL

Sample Solvent: Mobile phase A

Response Time: 0.12 sec

Flow Cell: 5.0 μL semi-micro

Data Rate: 14 Hz

LC System: Agilent 1100 Quaternary

